



INDICATIONS FOR USE

Aurstat™ Skin & Wound HydroGel is intended for use to relieve itch and pain from minor skin irritations, lacerations, abrasions and minor burns.

WARNING

For external use only. Not for injection. Do not use in the eyes. In case of deep or puncture wounds, animal bites or serious burns, consult a health care professional. Treatment of skin ulcers should be under the supervision of a health care professional. Stop use and consult a health care professional (1) if the condition persists, (2) signs of infection occur (e.g., increased pain, bleeding, wound drainage), (3) there is a change in wound color and/or odor, (4) the wound continues to get larger after the first few treatments, and (5) the wound does not begin to show signs of healing. If condition persists for longer than 7 days, consult a health care professional. Keep out of reach of children.

HANDLING

Non-irritating, non-cytotoxic and non-sensitizing to skin and eyes. No special handling precautions required.

STORAGE

Store with cap closed tightly after each use in its original container at room temperature away from direct sunlight and heat. Do not freeze.

DISPOSAL

No special disposal considerations required.

DIRECTIONS FOR USE

Cleanse the area and blot it dry. Apply a thin layer of Aurstat™ Skin and Wound HydroGel to the affected and immediately surrounding area. Use once or twice daily or as directed by physician. Cover with a bandage if necessary.

INGREDIENTS

Electrolyzed Water (96.524%), Sodium Magnesium Fluorosilicate (3.000%), Sodium Phosphate (0.400%), Sodium Chloride (0.066%), Hypochlorous Acid (0.008%), and Sodium Hypochlorite (0.002%)

Manufactured for:



Onset Dermatology
Cumberland, RI 02864
(888) 713-8154
www.onsetdermatology.com

Patent Protected

P/N: 2148-pf Rev. 2W 07-11